

BA-400

AUTOMATED CLINICAL CHEMISTRY ANALYZER



AUTOMATED CHEMISTRY ANALYZER

Smart Efficiency

BioSystems designs and develops efficient systems that implement the latest and best technical solutions.

BioSystems' BA400 is a clinical chemistry and turbidimetry analyzer designed to offer the best performance to laboratories looking towards achieving highest efficiency with optimal operative cost.

In combination with BioSystems original reagents and worldwide technical support coverage, the BA400 system defines the new generation of clinical analyzers.

Smart Autonomy

88 refrigerated positions with internal barcode reader. 135 positions for samples, controls and standards suitable for primary or paediatric tubes, 90 of them with barcode reading.

High capacity washing solution and waste containers, able to operate up to 8 hours of continuous working without refilling/voiding.

Automated water inlet and waste outlet with internal reservoirs, easy to adapt to any lab facilities.



Biosystems has developed for its BA400 analyzer an advanced and innovative optical system based on a battery of high power LED monochromatic sources with 8 working wavelengths that covers the most demanding methods of routine and special chemistry.

Solid-state light source with a split reference beam, with a working life up to 50.000 hours, to achieve optimal accuracy and performance.


Smart Performance

Self-controlled electronic subsystems through CAN bus optimize performance and reduce maintenance down-times. Sample dispensing system of high accuracy with level, collision and clot detection that automatically retreats to a protected home position during stops.

Low water consumption (less than 14 I/hour) with thermostated fluid washing station system to keep rotor temperature stable.

Low mechanical wear dispensing pumps with ceramic piston.

Independent powered cooling system for reagents (temperature between 5 and 8 oC, up to 35 oC room temperature).










Smart Solutions


High operating autonomy, through its high capacity for samples and reagents. Optical system with Biosystems' patented LED technology, with virtually no maintenance.

Low operating cost with optimized water and power consumption, minimum reaction volume and high pipetting precision.

Distributed electronics through CAN (Controller Area Network) bus system to increase robustness, simplify maintenance and reduce down times.

User friendly software, with intuitive graphical interface, real time monitoring of work-session and exhaustive quality control analysis (Westgard rules, Youden and Levy-Jennings charts, historical results database management).







Smart System

Original reagents specially designed and optimized for its use in the BA400 system, covering a complete panel of clinical chemistry and specific protein tests.

Worldwide technical assistance coverage with Biosystems' certified engineers. BioSystems SA, as a European manufacturer of its own reagents and analyzers, ensures proper functionality of all components under strict quality and safety standards for maximum performance and capabilities of their systems.





Technical Specifications

135

Yes

90

0.1 μL

Yes

Yes

1 uL

Yes

600 μL 120

± 0.2 °C ± 0.1 °C

UV methacrylate 37 ℃

7 tips (2 wash, 3 rinse, 2 dry)

CV < 1% to 0.1 A

CV < 0.1% to 2 A

Throughput ISE module Principles of analysis

400 test/h (without electrolytes) 320 test/h Colorimetry, turbidimetry.

ISE direct Serum, Plasma or Urine

(height up to 100 mm)

13.5 mm diameter

from 2 µL to 40 µL

From 1:2 to 1:200

20 mL, 60 mL

44 bottles of 20 mL)

From 90 µL to 450 µL From 10 µL to 300 µL

From 5 °C to 8 °C (at room temperature of 25 °C)

Na+, K+, Cl-. Li+ (optional) Serum: 100 μL / Urine: 200 μL

Diameter from 12 mm to 16 mm

Ceramic piston pump with low maintenance

Serum, Plasma, Whole Blood, Urine Cerebrospinal, Semen, Biological Liquid

88 (44 bottles of 20 mL or 60 mL +

Ceramic piston pump with low-maintenance

ISE Module (optional)

Method type Sample type Electrode type Sample volume

Sample handling

Capacity of sample rotor Barcode detector Number of samples with barcode Sample tube size

Pediatric well Type of syringe Sample type:

Pipetting volume Pipetting resolution Predilution ratio Level detection Clot detector Vertical collision detector

Reagent handling

Volume of reagent bottles Capacity of reagent rotor

Temperature range of refrigerator

Barcode detector Reagent volume R1 Reagent volume R2 Type of syringe Pipetting resolution Level detection Vertical collision detector Thermostated tip

Reactions rotor

Minimum reaction volume 1 Maximum reaction volume Number of wells Well material Temperature reaction rotor Accuracy of temperature Temperature stability

Cuvette washing system

Optical System

Light source Wavelengths Filters bandwidth Photometric range Internal resolution Detector Measurement precision (for 340 nm,

405 nm and 505 nm)

LED + Hard Coating filter 340 - 405 - 505 - 535 - 560 - 600 - 635 - 670 nm 10 nm + 2 nm-0.2 A to 3.5 A 0.0001 A Main Photodiode + reference photodiode

Environmental Requirement

Ambient temperature From 10 °C to 35 °C From 10 °C to 30 °C (With ISE module)

Relative humidity < 85% without condensation

< 2 500 m Altitude

Dimensions and weight

Dimensions (width, depth and height) 1 200 mm x 720 mm x 1 258 mm Weight

Electrical Requirements

115 V to 230 V Main voltages Main frequencys 50 Hz or 60 Hz Electric power 500 VA

Fluidic Requirements

Water inlet External tank or mains water supply

Water type Purified Type II Water consumption < 14 L/h Bottle of high concentration waste 5 L Bottle of washing solution 5 L

Minimum Computer Requirements

Windows® 7, 8.1 or 10 Operating system

CPU Equivalent to Intel Core i3 @3.10 GHz or higger

RAM 4 Gbytes 40 GB or higher Hard Disk

DVD Yes

Minimum resolution 1024 x 768 Monitor minimum resolution

Connector of serial channel

Connectivity to LIS HL7 and ASTM protocols

(Laboratory Information Systems)

Directives and Standards Compliance

EC Directive 98/79/EC IVD



Biochemistry

Code	Test	Presentation	
		R	R2
21550	α-Amylase-Direct	8x20 mL	
21534	α-Amylase-EPS	2x60 mL	2x15 mL
21799	α-Amylase Pancreatic	2x60 mL	2x15 mL
2153	Alanine Aminotransferase (ALT/GPT)	8x60 mL	8x15 mL
21547	Albumin	10x60 mL	
21592	Alkaline Phosphatase (ALP)-AMP	4x60 mL	4x15 mL
21590	Alkaline Phosphatase (ALP)-DEA	4x60 mL	4x15 mL
21531	Aspartate Aminotransferase (AST/GOT)	8x60 mL	8x15 mL
21798	Bilirubin (Direct)	4x60 mL	3x20 mL
21510	Bilirubin (Total)	8x60 mL	8x15 mL
21570	Calcium-Arsenazo	10x60 mL	
21558	Carbon Dixide	2x60 mL	
21505	Cholesterol	10x60 mL	
21557	Cholesterol HDL Direct	2x60 mL	2x20 mL
21585	Cholesterol LDL Direct	2x60 mL	2x20 mL
21588	Cholinesterase (CHE)	2x60 mL	2x15 mL
21790	Creatine Kinase (CK)	2x60 mL	2x15 mL
21792	Creatine Kinase-MB (CK-MB)	2x60 m	2x15 mL
21502	Creatinine	5x60 mL	5x60 mL
21520	γ-Glutamyltransferase (γ-GT)	4x60 mL	4x15 mL
21503	Glucose	10x60 mL	
21509	Iron-Ferrozine	4x60 mL	4x15 mL
21580	Lactate Dehydrogenase (LDH)	8x60 mL	8x15 mL
21586	Lactate Dehydrogenase (LDH)-IFCC	8x60 mL	8x15 mL
21793	Lipase	2x50 mL	1x20 mL
21797	Magnesium	2x60 mL	2x15 mL
21508	Phosphorus	4x60 mL	2x50 mL
21500	Protein (Total)	10x60 mL	
21501	Protein (Urine+CSF)	4x60 mL	
21528	Triglycerides	10x60 mL	
21516	Urea/BUN-UV	8x60 mL	8x15 mL
21521	Uric Acid	10x60 mL	



Turbidimetry

Code	Test	Present	Presentation	
		R	R2	
22324	Albumin (Microalbuminuria)	4x60 mL	4x15 mL	
22923	Anti-Streptolysin O (ASO)	2x60 mL	2x15 mL	
22936	Antithrombin III	2x60 mL	2x15 mL	
22928	α-1-Acid Glycoprotein	2x60 mL		
22491	α-1-Microglobulin	2x60 mL	2x15 mL	
22095	Apolipoprotein A-I (APO A-I)	2x60 mL	2x15 mL	
22098	Apolipoprotein B (APO B)	2x60 mL	2x15 mL	
22925	b2-Microglobulin	2x60 mL	2x15 mL	
22084	Complement Component C3	2x60 mL		
22085	Complement Component C4	2x60 mL		
22921	C-Reactive Protein (CRP)	4x60 mL	4x15 mL	
22927	C-Reactive Protein-hs (CRP-hs)	2x60 mL	2x15 mL	
22044	Hemoglobin A1C-turbi	75 mL		
22934	Ferritin	2x50 mL	1x50 mL	
22082	Immunoglobulin A (lg A)	2x60 mL		
22081	Immunoglobulin G (lg G)	2x60 mL		
22083	Immunoglobulin M (Ig M)	2x60 mL		
22929	Prealbumin	2x60 mL		
22922	Rheumatoid Factors (RF)	4x60 mL	4x15 mL	
22091	Transferrin	2x60 mL		

